

Olerup SSP[®] DPA1

Product number:	101.331-24/06 – including <i>Taq</i> pol. 101.331-24u/06u – without <i>Taq</i> pol.
Lot number:	73G
Expiry date:	2012-January-01
Number of tests:	24 test – Product No. 101.331-24/-24u 6 tests – Product No. 101.331-06/-06u
Number of Wells per test:	16

CHANGES COMPARED TO THE PREVIOUS DPA1 LOT (20F):

The DPA1 primer set is unchanged.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

DPA1*0103 to 0110, DPA1*0201 to 0204, DPA1*0301 to 0303 and DPA1*0401, i.e. all the currently recognized DPA1 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2009-October-19, release 2.27.0.

RESOLUTION IN DPA1 HOMO- AND HETEROZYGOTES:

Excellent. All frequently occurring genotypes give rise to unique amplification patterns.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.