

## **Olerup SSP<sup>®</sup> DPA1**

Product number:	101.331-24/06 – including <i>Taq</i> pol. 101.331-24u/06u – without <i>Taq</i> pol.
Lot number:	53K
Expiry date:	2012-September-01
Number of tests:	24 test – Product No. 101.331-24/24u 6 tests – Product No. 101.331-06/06u
Number of wells per test:	16

### **CHANGES COMPARED TO THE PREVIOUS DPA1 LOT (73G):**

The DPA1 primer set is unchanged.

**THE NUMBER OF WELLS** is unchanged.

### **ALLELE COVERAGE:**

DPA1\*01:03 to DPA1\*01:10, DPA1\*02:01 to DPA1\*02:04, DPA1\*03:01 to DPA1\*03:03 and DPA1\*04:01, i.e. all the currently recognized DPA1 alleles, give rise to unique amplification patterns; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2010-July-16, release 3.1.0.

### **RESOLUTION IN DPA1 HOMO- AND HETEROZYGOTES:**

Excellent. All frequently occurring genotypes give rise to unique amplification patterns.

### **MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No comments received.