

Olerup SSP[®] DPB1

Product number: 101.311-24/04 – including *Taq* pol.
101.311-24u/04u – without *Taq* pol.
Lot number: 09M
Expiry date: 2013-November-01
Number of tests: 24 tests – Product No. 101.311-24
4 tests – Product No. 101.311-04
Number of tubes per test: 48

CHANGES COMPARED TO THE PREVIOUS DPB1 LOT (32K):

Well	5'-primer	3'-primer	rationale
16	Added	Added	Primer pair added for the DPB1*127:01 allele.
23	Added	Added	Primer pair added for the DPB1*127:01 allele.
32	Added	-	Primer added for the DPB1*130:01 allele.
34	-	-	Exchanged positive control primer pair.
37	-	-	Exchanged positive control primer pair.
42	-	Added	Primer added for the DPB1*128:01 allele.
46	-	-	Exchanged positive control primer pair.
47	-	-	Exchanged positive control primer pair.
48	-	-	Exchanged positive control primer pair.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

DPB1*01:01 to DPB1*132:01, i.e. all the currently recognized DPB1 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2011-January-14, release 3.3.0.

RESOLUTION IN DPB1 HOMO- AND HETEROZYGOTES:

The DPB1*01:01,01:01, DPB1*01:01,89:01 and DPB1*01:01,90:01 give rise to identical amplification patterns.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.