

Olerup SSP[®] DPB1

Product number: 101.311-24/04 – including *Taq* pol.
101.311-24u/04u – without *Taq* pol.
Lot number: 09N
Expiry date: 2014-August-01
Number of tests: 24 tests – Product No. 101.311-24/24u
4 tests – Product No. 101.311-04/04u
Number of tubes per test: 48

CHANGES COMPARED TO THE PREVIOUS DPB1 LOT (09M):

Well	5'-primer	3'-primer	rationale
3	Modified	-	Improved specificity of primer pair.
41	-	Added	Primer added for the DPB1*134:01 allele.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

DPB1*01:01 to DPB1*136:01, i.e. all the currently recognized DPB1 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2012-January-12, release 3.7.0.

RESOLUTION IN DPB1 HOMO- AND HETEROZYGOTES:

The DPB1*01:01,01:01, DPB1*01:01,89:01 and DPB1*01:01,90:01 give rise to identical amplification patterns.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.