

## **Olerup SSP<sup>®</sup> DPB1**

Product number:	101.311-24/04 – including <i>Taq</i> pol. 101.311-24u/04u – without <i>Taq</i> pol.
Lot number:	12R
Expiry date:	2015-March-01
Number of tests:	24 tests – Product No. 101.311-24/24u 4 tests – Product No. 101.311-04/04u
Number of tubes per test:	48

### **CHANGES COMPARED TO THE PREVIOUS DPB1 LOT (09N):**

The DPB1 primer set is unchanged compared to the previous lot.

**THE NUMBER OF WELLS** is unchanged.

### **ALLELE COVERAGE:**

DPB1\*01:01 to DPB1\*139:01, i.e. all the currently recognized DPB1 alleles, give rise to unique amplification patterns; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2012-July-12, release 3.9.0.

### **RESOLUTION IN DPB1 HOMO- AND HETEROZYGOTES:**

The DPB1\*01:01,01:01, DPB1\*01:01,89:01 and DPB1\*01:01,90:01 give rise to identical amplification patterns.

### **MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No comments received.