

Lot No.: **30K**

Olerup SSP[®] DQ low resolution

Product number: 101.201-48/12 - including *Taq* pol.
101.201-48u/12u – without *Taq* pol.
Lot number: 30K
Expiry date: 2012-July-01
Number of tests: 48 tests – Product No. 101.201-48
12 tests – Product No. 101.201-12
Number of wells per test: 8

Changes compared to the previous DQ low resolution Lot (91F):

Well	5'-primer	3'-primer	rationale
4	Removed, added	Removed, added	Primer pair from well 5.
5	Exchanged	Exchanged	Improved resolution of the DQ7, DQ8 and DQ9 specificities.
6	-	Exchanged	Improved resolution of the DQ7, DQ8 and DQ9 specificities.
7	Added	-	Primer added for the DQB1*03:25 allele.
8	Added	-	Primer added for the DQB1*04:03:02 allele.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

DQB1*05:01:01 to 05:05, DQB1*06:01:01 to 06:40, DQB1*02:01:01 to 02:05, DQB1*03:01:01 to 03:26 and DQB1*04*01 to 04:04, i.e. all the currently recognized DQB1 alleles will be grouped into their corresponding serological specificities; www.ebi.ac.uk/imgt/hla, 2010-April-01, release 3.0.0.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.