

101.213-24 – including *Taq* polymerase101.213-24u – without *Taq* polymeraseLot No.: **66E****Olerup SSP<sup>®</sup> DQB1\*02**

**Product number:** 101.213-24 – including *Taq* polymerase  
101.213-24u – without *Taq* polymerase  
**Lot number:** 66E  
**Expiry date:** 2010-June-01  
**Number of tests:** 24  
**Number of Wells per test:** 6

**CHANGES COMPARED TO THE PREVIOUS DQB1\*02 LOT (X56):**

The DQB1\*02 primer set is unchanged.

THE NUMBER OF WELLS is unchanged.

**ALLELE COVERAGE:**DQB1\*0201 to DQB1\*0205, i.e. all the currently recognized DQB1\*02 alleles, give rise to unique amplification patterns; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2008-April-08, release 2.21.0.**RESOLUTION IN DQB1\*02 HOMO- AND HETEROZYGOTES:**

Very good.

**INFLUENCE ON THE INTERPRETATION OF DQB1\*02 SUBTYPINGS BY NON-DQB1\*02 ALLELES:**

None frequently occurring.

**MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No suggestions received.