

101.213-24 – including *Taq* polymerase101.213-24u – without *Taq* polymeraseLot No.: **08G****Olerup SSP® DQB1*02**

Product number: 101.213-24 – including *Taq* polymerase
101.213-24u – without *Taq* polymerase
Lot number: 08G
Expiry date: 2011-August-01
Number of tests: 24
Number of Wells per test: 6

CHANGES COMPARED TO THE PREVIOUS DQB1*02 LOT (66E):

The DQB1*02 primer set is unchanged.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:DQB1*0201 to DQB1*0205, i.e. all the currently recognized DQB1*02 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2009-July-17, release 2.26.0.**RESOLUTION IN DQB1*02 HOMO- AND HETEROZYGOTES:**

Very good.

INFLUENCE ON THE INTERPRETATION OF DQB1*02 SUBTYPINGS BY NON-DQB1*02 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.