

101.213-24 – including *Taq* polymerase101.213-24u – without *Taq* polymeraseLot No.: **77K****Olerup SSP® DQB1*02**

Product number:	101.213-24 – including <i>Taq</i> polymerase 101.213-24u – without <i>Taq</i> polymerase
Lot number:	77K
Expiry date:	2013-July-01
Number of tests:	24
Number of Wells per test:	6

CHANGES COMPARED TO THE PREVIOUS DQB1*02 LOT (08G):

The DQB1*02 primer set is unchanged.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:DQB1*02:01 to DQB1*02:05, i.e. all the currently recognized DQB1*02 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2010-October-15, release 3.2.0.**RESOLUTION IN DQB1*02 HOMO- AND HETEROZYGOTES:**

Very good.

INFLUENCE ON THE INTERPRETATION OF DQB1*02 SUBTYPINGS BY NON-DQB1*02 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.