

DQB1*03
101.214-24/06 – licensed for PCR
101.214-24u/06u – not licensed for PCR
Lot No.: **Y07**

1

Release Note - *Olerup* SSP™ DQB1*03 – Lot No. **Y07**

Product number: 101.214-24/06 – licensed for PCR
101.214-24u/06u – not licensed for PCR
Lot number: Y07
Expiry date: 2009-July-01
Number of tests: 24 test – Product No. 101.214-24
6 tests – Product No. 101.214-06
Number of tubes per test: 16

CHANGES COMPARED TO THE PREVIOUS DQB1*03 Lot (v55):

Tube	5'-primer	3'-primer	Rationale
4	Added	Added	Primer pair added for the DQB1*0319 allele.
6	Added	-	Primer added for the DQB1*0320 allele.
14	-	Added	Primer added for the DQB1*030103 allele.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

DQB1*0301 to DQB1*0320, i.e. all the currently recognized DQB1*03 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2007-July-09, release 2.7.

RESOLUTION IN DQB1*03 HOMOZYGOTES:

The DQB1*0301,0302 and DQB1*0302,0319 genotypes give rise to the same amplification patterns.

The DQB1*0301,0303, DQB1*0301,0319, DQB1*0303,0319 and DQB1*0319,0319 genotypes give rise to the same amplification patterns.

INFLUENCE ON THE INTERPRETATION OF DQB1*03 SUBTYPINGS BY NON-DQB1*03 ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions have been received.