

101.214-24/06 – including *Taq* polymerase
101.214-24u/06u – without *Taq* polymerase

Lot No.: **22G**

Olerup SSP® DQB1*03

Product number: 101.214-24/06 – including *Taq* pol.
101.214-24u/06u – without *Taq* pol.
Lot number: 22G
Expiry date: 2011-October-01
Number of tests: 24 test – Product No. 101.214-24
6 tests – Product No. 101.214-06
Number of wells per test: 23

CHANGES COMPARED TO THE PREVIOUS DQB1*03 LOT (15F):

Well	5'-primer	3'-primer	rationale
22	New	New	New primer pair for the DQB1*0323 allele.
23	New	New	New primer pair for the DQB1*0324 allele.

THE NUMBER OF WELLS has been increased from 21 to 23.

ALLELE COVERAGE:

DQB1*0301 to DQB1*0325, i.e. all the currently recognized DQB1*03 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2009-July-17, release 2.26.0.

RESOLUTION IN DQB1*03 HOMOZYGOTES:

The DQB1*0301,0301, DQB1*0301,0302, DQB1*0302,0302, DQB1*0302,0303 and DQB1*0303,0303 genotypes all give rise to unique amplification patterns.

INFLUENCE ON THE INTERPRETATION OF DQB1*03 SUBTYPINGS BY NON-DQB1*03 ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions have been received.