

DQB1*04
101.215-12 – licensed for PCR
101.215-12u – not licensed for PCR
Lot No.: Y50

1

Release Note - *Olerup SSP*TM DQB1*04 – Lot No. Y50

Product number:	101.215-12 – licensed for PCR 101.215-12u – <u>not</u> licensed for PCR
Lot number:	Y50
Expiry date:	2009-September-01
Number of tests:	12
Number of tubes per test:	3

CHANGES COMPARED TO THE PREVIOUS DQB1*04 LOT (V80):

The DQB1*04 primer set is unchanged.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

DQB1*0401 to DQB1*0402, i.e. all the currently recognized DQB1*04 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2007-October-05, release 2.19.0.

RESOLUTION IN DQB1*04 HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF DQB1*04 SUBTYPINGS BY NON-DQB1*04 ALLELES:

None.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No modifications made.