

Olerup SSP[®] DQB1*04

Product number: 101.215-12 – including *Taq* polymerase
101.215-12u – without *Taq* polymerase
Lot number: 67E
Expiry date: 2010-June-01
Number of tests: 12
Number of wells per test: 4

CHANGES COMPARED TO THE PREVIOUS DQB1*04 LOT (Y50):

Well	5'-primer	3'-primer	rationale
1	Exchanged	-	Primer exchanged to amplify all DQB1*04 alleles.
4	New	New	New primer pair for the DQB1*0403 allele.

THE NUMBER OF WELLS has been increased from 3 to 4.

ALLELE COVERAGE:

DQB1*0401 to DQB1*0403, i.e. all the currently recognized DQB1*04 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2008-April-08, release 2.21.0.

RESOLUTION IN DQB1*04 HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF DQB1*04 SUBTYPINGS BY NON-DQB1*04 ALLELES:

None.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No modifications made.