

101.215-12 – including *Taq* polymerase101.215-12u – without *Taq* polymeraseLot No.: **14G****Olerup SSP® DQB1*04**

Product number: 101.215-12 – including *Taq* polymerase
101.215-12u – without *Taq* polymerase
Lot number: 14G
Expiry date: 2011-July-01
Number of tests: 12
Number of wells per test: 4

CHANGES COMPARED TO THE PREVIOUS DQB1*04 LOT (67E):

The DQB1*04 primer set is unchanged.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:DQB1*0401 to DQB1*0403, i.e. all the currently recognized DQB1*04 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2009-May-11, release 2.25.2.**RESOLUTION IN DQB1*04 HOMO- AND HETEROZYGOTES:**

Excellent.

INFLUENCE ON THE INTERPRETATION OF DQB1*04 SUBTYPINGS BY NON-DQB1*04 ALLELES:

None.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No modifications made.