

## Olerup SSP® DQB1\*04

**Product number:** 101.215-12 – including *Taq* polymerase  
101.215-12u – without *Taq* polymerase  
**Lot number:** 54K  
**Expiry date:** 2012-September-01  
**Number of tests:** 12  
**Number of wells per test:** 5

### CHANGES COMPARED TO THE PREVIOUS DQB1\*04 LOT (14G):

Well	5'-primer	3'-primer	rationale
1	Added	-	Primer added for the DQB1*04:03:02 allele.
5	New	New	New primer pair for the DQB1*04:04 and DQB1*04:05 alleles.

**THE NUMBER OF WELLS** is increased from 4 to 5.

#### ALLELE COVERAGE:

DQB1\*04:01 to DQB1\*04:05, i.e. all the currently recognized DQB1\*04 alleles, give rise to unique amplification patterns<sup>1</sup>; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2010-July-16, release 3.1.0.

<sup>1</sup>The DQB1\*04 subtyping kit cannot distinguish the DQB1\*04:01:01 and DQB1\*04:01:02 alleles.

#### RESOLUTION IN DQB1\*04 HOMO- AND HETEROZYGOTES:

Excellent.

#### INFLUENCE ON THE INTERPRETATION OF DQB1\*04 SUBTYPINGS BY NON-DQB1\*04 ALLELES:

None.

#### MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No modifications made.