

Olerup SSP® DQB1*04

Product number:	101.215-12 – including <i>Taq</i> polymerase 101.215-12u – without <i>Taq</i> polymerase
Lot number:	62R
Expiry date:	2015-May-01
Number of tests:	12
Number of wells per test:	8

CHANGES COMPARED TO THE PREVIOUS DQB1*04 LOT (87M):

The DQB1*04 primer set is unchanged compared to the previous Olerup SSP® DQB1*04 lot.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

DQB1*04:01 to DQB1*04:08, i.e. all the currently recognized DQB1*04 alleles, give rise to unique amplification patterns¹; www.ebi.ac.uk/imgt/hla, 2012-July-12, release 3.9.0.

¹The DQB1*04 subtyping kit cannot distinguish the silent mutations in the DQB1*04:01:01 and DQB1*04:01:02 alleles or the DQB1*04:02:01 and 04:02:02 alleles.

RESOLUTION IN DQB1*04 HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF DQB1*04 SUBTYPINGS BY NON-DQB1*04 ALLELES:

None.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No modifications made.