

DQB1\*05  
101.211-24 – including *Taq* polymerase  
101.211-24u – without *Taq* polymerase  
Lot No.: **09E**

Release Note

Page 1 of 1

## **Olerup SSP<sup>®</sup> DQB1\*05**

**Product number:** 101.211-24 – including *Taq* polymerase  
101.211-24u – without *Taq* polymerase  
**Lot number:** 09E  
**Expiry date:** 2010-February-01  
**Number of tests:** 24  
**Number of Wells per test:** 6

### **CHANGES COMPARED TO THE PREVIOUS DQB1\*05 LOT (X15):**

The DQB1\*05 primer set is unchanged.

**THE NUMBER OF WELLS** is unchanged.

### **ALLELE COVERAGE:**

DQB1\*0501 to DQB1\*0505, i.e. all the currently recognized DQB1\*05 alleles, give rise to unique amplification patterns; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2008-January-11, release 2.20.0.

### **RESOLUTION IN DQB1\*05 HOMO- AND HETEROZYGOTES:**

Very good.

### **INFLUENCE ON THE INTERPRETATION OF DQB1\*05 SUBTYPINGS BY NON-DQB1\*05 ALLELES:**

None.

### **MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No suggestions received.