

## Olerup SSP<sup>®</sup> DQB1\*05

Product number: 101.211-24 – including *Taq* polymerase  
101.211-24u – without *Taq* polymerase  
Lot number: 66R  
Expiry date: 2015-June -01  
Number of tests: 24  
Number of wells per test: 15

### CHANGES COMPARED TO THE PREVIOUS DQB1\*05 LOT (83N):

Well	5'-primer	3'-primer	rationale
13	New	New	New primer pair for the DQB1*05:15 allele.
14	New	New	New primer pair for the DQB1*05:16 allele.
15	New	New	New primer pair for the DQB1*05:17 allele.

THE NUMBER OF WELLS has been increased from 12 to 15.

### ALLELE COVERAGE:

DQB1\*05:01 to DQB1\*05:18, i.e. all the currently recognized DQB1\*05 alleles, give rise to unique amplification patterns; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2012-October-17, release 3.10.0.

### RESOLUTION IN DQB1\*05 HOMO- AND HETEROZYGOTES:

Very good.

### INFLUENCE ON THE INTERPRETATION OF DQB1\*05 SUBTYPINGS BY NON-DQB1\*05 ALLELES:

None.

### MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.