

Olerup SSP[®] DQB1*05

Product number:	101.211-24 – including <i>Taq</i> polymerase 101.211-24u – without <i>Taq</i> polymerase
Lot number:	4S6
Expiry date:	2027-09-01
Number of tests:	24
Number of wells per test:	31+1

Changes compared to the previous DQB1*05 Lot (0R1):

The DQB1*05 resolution primer set is unchanged compared to the previous *Olerup SSP[®] DQB1*05* lot (Lot No. 0R1).

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

DQB1*05:01 to DQB1*05:323, i.e. all the currently recognized DQB1*05 alleles, will be amplified by the primers in the DQB1*05 subtyping kit¹, www.ebi.ac.uk/imgt/hla, 2023-July-12, release 3.53.0.

The DQB1*05 kit enables separation of the confirmed DQB1*05 alleles as listed in the IMGT/HLA database 3.26.0. An HLA allele is listed as confirmed by IMGT/HLA if it has been sequenced by more than a single laboratory or from multiple sources.

The DQB1*05 also enables identification of many null and alternatively expressed alleles.

The following DQB1*05 alleles can be distinguished by the different sizes of the HLA-specific PCR product:

Alleles	Primer mix	Alleles	Primer mix
DQB1*05:12, 05:242	10	DQB1*05:27, 05:29	8
DQB1*05:13, 05:201, 05:259	11	DQB1*05:35, 05:117	11
DQB1*05:16, 05:101	14	DQB1*05:45, 05:169, 05:182	11
DQB1*05:20:01, 05:159	16	DQB1*05:110N, 05:156	27

The following alleles will give rise to identical amplification patterns with the DQB1*05 subtyping kit. These alleles can be distinguished by the DQ low resolution kit and/or the DQB1*06 subtyping kit:

Alleles

DQB1*05:300, DQB1*06:389

¹Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.



DQB1*05

101.211-24 – including *Taq* polymerase

101.211-24u – without *Taq* polymerase

Lot No.: **4S6**

RESOLUTION IN DQB1*05 HOMO- AND HETEROZYGOTES:

Very good.

INFLUENCE ON THE INTERPRETATION OF DQB1*05 SUBTYPINGS BY NON-DQB1*05 ALLELES:

None.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.

