

101.212-24 – including *Taq* polymerase101.212-24u – without *Taq* polymeraseLot No.: **61E****Olerup SSP<sup>®</sup> DQB1\*06**

**Product number:** 101.212-24 – including *Taq* polymerase  
 101.212-24u – without *Taq* polymerase  
**Lot number:** 61E  
**Expiry date:** 2010-May-01  
**Number of tests:** 24  
**Number of Wells per test:** 28

**CHANGES COMPARED TO THE PREVIOUS DQB1\*06 LOT.**

Well	5'-primer	3'-primer	rationale
10	Removed	-	Primer pair moved to well 26.
23	-	Removed	Primer pair moved to well 25.
25	New	New	Primer pair moved from well 23.
26	New	New	Primer pair moved from well 10
27	New	New	Primer pair added for the DQB1*0633 allele
28	New	New	Primer pair added for the DQB1*0634 allele

**THE NUMBER OF WELLS** is increased from 24 to 28.

**ALLELE COVERAGE:**

**DQB1\*0601 to DQB1\*0634** i.e. all the currently recognized DQB1\*06 alleles, give rise to unique amplification patterns; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2008-April-08, release 2.21.0.

**RESOLUTION IN DQB1\*06 HOMO- AND HETEROZYGOTES:**

Very good.

**INFLUENCE ON THE INTERPRETATION OF DQB1\*06 SUBTYPINGS BY NON-DQB1\*06 ALLELES:**

None.

**MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No suggestions received.