

Olerup SSP[®] DQA1*02,05;DQB1*02,0302

Product number: 101.903-24– including *Taq* pol.
101.903-24u– without *Taq* pol.
Lot number: 33G
Expiry date: 2011-September-01
Number of tests: 24 tests – Product No. 101.903-24
Number of wells per test: 16

CHANGES COMPARED TO THE PREVIOUS DQA1*02,05;DQB1*02,0302 LOT (07F):

Well	5'-primer	3'-primer	rationale
11	Added	-	Primers added for the DQB1*0323 and 0325 alleles.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

All the currently recognized DQA1*02, DQA1*05 and DQB1*02 and DQB1*0302 alleles give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2009-July-17, release 2.26.0.

RESOLUTION IN HOMOZYGOTES:

Good.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.