

Olerup SSP[®] DQA1*02,05;DQB1*02,03:02

Product number: 101.903-24– including *Taq* pol.
101.903-24u– without *Taq* pol.
Lot number: 50K
Expiry date: 2012-September-01
Number of tests: 24 tests – Product No. 101.903-24
Number of wells per test: 15+1

CHANGES COMPARED TO THE PREVIOUS DQA1*02,05;DQB1*02,03:02 Lot (33G):

Well	5'-primer	3'-primer	rationale
6	-	Exchanged	3'-primer exchanged for improved specificity.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

All the currently recognized DQA1*02, DQA1*05 and DQB1*02 and DQB1*03:02 alleles give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2010-July-16, release 3.1.0.

RESOLUTION IN HOMOZYGOTES:

Good.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.