

## Olerup SSP<sup>®</sup> DQA1\*02,05;DQB1\*02,03:02

**Product number:** 101.903-24– including *Taq* pol.  
101.903-24u– without *Taq* pol.  
**Lot number:** 34M  
**Expiry date:** 2013-November-01  
**Number of tests:** 24 tests – Product No. 101.903-24  
**Number of wells per test:** 15+1

### CHANGES COMPARED TO THE PREVIOUS DQA1\*02,05;DQB1\*02,03:02 Lot (50K):

Well	5'-primer	3'-primer	rationale
8	Added	-	5'-primer added for the DQB1*03:37 allele.
9	Added	Added	Primer pair added for the DQB1*03:32 allele.

**THE NUMBER OF WELLS** is unchanged.

#### ALLELE COVERAGE:

All the currently recognized DQA1\*02, DQA1\*05 and DQB1\*02 and DQB1\*03:02 alleles give rise to unique amplification patterns; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2011-April-08, release 3.4.0.

#### RESOLUTION IN HOMOZYGOTES:

Good.

#### MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.