

## **Olerup SSP<sup>®</sup> DQA1\*02,05;DQB1\*02,03:02**

**Product number:** 101.903-24– including *Taq* pol.  
101.903-24u– without *Taq* pol.  
**Lot number:** 59N  
**Expiry date:** 2014-September-01  
**Number of tests:** 24 tests – Product No. 101.903-24  
**Number of wells per test:** 15+1

### **CHANGES COMPARED TO THE PREVIOUS DQA1\*02,05;DQB1\*02,03:02 Lot (34M):**

<b>Well</b>	<b>5'-primer</b>	<b>3'-primer</b>	<b>rationale</b>
6	-	Exchanged	Exchanged 3'-primer for improved resolution of the DQB1*03:02 and DQB1*03:04 alleles.

**THE NUMBER OF WELLS** is unchanged.

#### **ALLELE COVERAGE:**

All the currently recognized DQA1\*02, DQA1\*05 and DQB1\*02 and DQB1\*03:02 alleles give rise to unique amplification patterns; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2012-January-12, release 3.7.0.

#### **RESOLUTION IN HOMOZYGOTES:**

Good.

#### **MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No comments received.