

Olerup SSP® DQA1*02,05;DQB1*02,0302

Product number: 101.903-24– including *Taq* pol.
101.903-24u– without *Taq* pol.
Lot number: 07F
Expiry date: 2010-October-01
Number of tests: 24 tests – Product No. 101.903-24
Number of wells per test: 16

CHANGES COMPARED TO THE PREVIOUS DQA1*02,05;DQB1*02,0302 Lot (05E):

Well	5'-primer	3'-primer	rationale
12	Moved	Moved	Primer pair from well 13
13	Moved	Moved	Primer pair from well 14
14	Moved	Moved	Primer pair from well 15
15	Moved	Moved	Primer pair from well 16
16	Moved	Moved	Negative control moved from well 12.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

All the currently recognized DQA1*02, DQA1*05 and DQB1*02 and DQB1*0302 alleles give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2008-July-11, release 2.22.0.

RESOLUTION IN HOMOZYGOTES:

Good.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.