

DQA1  
101.231-24/04 – licensed for PCR  
101.231-24u/04u – not licensed for PCR  
Lot No.: **Y46**

1

## Release Note - *Olerup* SSP™ DQA1 – Lot No. **Y46**

**Product number:** 101.231-24/04 – licensed for PCR  
101.231-24u/04u – not licensed for PCR  
**Lot number:** Y46  
**Expiry date:** 2009-October-01  
**Number of tests:** 24 tests – Product No. 101.231-24  
4 tests – Product No. 101.231-04  
**Number of tubes per test:** 32

### CHANGES COMPARED TO THE PREVIOUS DQA1 LOT (X27):

Tube	5'-primer	3'-primer	rationale
31	Exchanged	-	Longer specific PCR product.

**THE NUMBER OF WELLS** is unchanged.

### ALLELE COVERAGE:

All the currently recognized DQA1 alleles, i.e. DQA1\*0101 to 0107, DQA1\*0201, DQA1\*0301 to 0303, DQA1\*0401 to DQA1\*0404, DQA1\*0501 to 0509 and DQA1\*0601 to DQA1\*0602, give rise to unique amplification patterns; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2007-October-05, release 2.19.0.

### RESOLUTION IN DQA1 HOMOZYGOTES:

Good.

### MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.