

Lot No.: 12F

Olerup SSP[®] DQB1 high resolution for frequent alleles

Product number: 101.221-12 – including *Taq* polymerase
101.221-12u – without *Taq* polymerase
Lot number: 12F
Expiry date: 2010-October-01
Number of tests: 12 tests
Number of wells per test: 47 + 1

CHANGES COMPARED TO THE PREVIOUS DQB1 HIGH RESOLUTION LOT (43E):

Well	5'-primer	3'-primer	rationale
4	Exchanged	-	Primer exchanged to decrease the tendency to form primer dimers.
9	Exchanged	-	Primer added for increased resolution.
37	Added	Added	Primer pair added for increased resolution
41	Added	Added	Primer pair added for the DQB1*0322 allele.
42	-	Added	Primer added for the DQB1*0403 allele.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

All the currently recognized DQB1 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2008-July-11, release 2.22.0.

RESOLUTION IN DQB1 HOMOZYGOTES AND HETEROZYGOTES:

Very good.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.