

101.901-24 – including *Taq* polymerase101.901-24u – without *Taq* polymeraseLot No.: **17G****Olerup SSP® DQB1*0602,DQA1*0102 - SSP**

Product number: 101.901-24 – including *Taq* polymerase
Product number: 101.901-24u – without *Taq* polymerase
Lot number: 17G
Expiry date: 2011-August-01
Number of tests: 24
Number of wells per test: 8

CHANGES COMPARED TO THE PREVIOUS DQB1*0602, DQA1*0102 Lot (37E):

Well	5'-primer	3'-primer	rationale
5	New	New	New primer pair for the DQA1*0106 allele.
6	New	New	New primer pair for the DQB1*0616 and *0619 alleles.
7	New	New	New primer pair for the DQB1*0620 and *0633 alleles.
8	New	New	New primer pair for the DQB1*0624 allele.

THE NUMBER OF WELLS has been increased from 4 to 8.

ALLELE COVERAGE:

DQB1*0601 to DQB1*0635 and DQA1*0101 to DQA1*0107, i.e. all the currently recognized DQB1*06 and DQA1*01 alleles, have been considered in the DQB1*0602,DQA1*0102 primer set; www.ebi.ac.uk/imgt/hla, 2009-July-17, release 2.26.0.

RESOLUTION IN DQB1*06 AND DQA1*01 HOMO- AND HETEROZYGOTES:

Very good.

INFLUENCE ON THE INTERPRETATION OF DQB1*0602, DQA1*0102 SUBTYPINGS BY OTHER ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No modifications made.