

101.901-24 – including *Taq* polymerase101.901-24u – without *Taq* polymeraseLot No.: **64M****Olerup SSP® DQB1*06:02,DQA1*01:02 - SSP**

Product number:	101.901-24 – including <i>Taq</i> polymerase
Product number:	101.901-24u – without <i>Taq</i> polymerase
Lot number:	64M
Expiry date:	2014-May-01
Number of tests:	24
Number of wells per test:	8

CHANGES COMPARED TO THE PREVIOUS DQB1*06:02, DQA1*01:02 Lot (96K):

The DQB1*06:02,DQA1*01:02 primer set is unchanged compared to the previous lot.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

DQB1*06:01 to DQB1*06:44 and DQA1*01:01 to DQA1*01:07, i.e. all the currently recognized DQB1*06 and DQA1*01 alleles, have been considered in the DQB1*06:02,DQA1*01:02 primer set; www.ebi.ac.uk/imgt/hla, 2011-July-14, release 3.5.0.

RESOLUTION IN DQB1*06 AND DQA1*01 HOMO- AND HETEROZYGOTES:

Very good.

INFLUENCE ON THE INTERPRETATION OF DQB1*06:02, DQA1*01:02 SUBTYPINGS BY OTHER ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No modifications made.