

Lot No.: **96K**

Olerup SSP[®] DQB1*06:02,DQA1*01:02 - SSP

Product number: 101.901-24 – including *Taq* polymerase
Product number: 101.901-24u – without *Taq* polymerase
Lot number: 96K
Expiry date: 2013-July-01
Number of tests: 24
Number of wells per test: 8

CHANGES COMPARED TO THE PREVIOUS DQB1*06:02, DQA1*01:02 Lot (17G):

Well	5'-primer	3'-primer	rationale
7	Added	-	Primers added for the DQB1*06:23 and *06:37 alleles.
8	-	Added	Primers added for the DQB1*06:14:01-06:14:02 and *06:15 alleles.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

DQB1*06:01 to DQB1*06:39 and DQA1*01:01 to DQA1*01:07, i.e. all the currently recognized DQB1*06 and DQA1*01 alleles, have been considered in the DQB1*06:02,DQA1*01:02 primer set; www.ebi.ac.uk/imgt/hla, 2010-October-15, release 3.2.0.

RESOLUTION IN DQB1*06 AND DQA1*01 HOMO- AND HETEROZYGOTES:

Very good.

INFLUENCE ON THE INTERPRETATION OF DQB1*06:02, DQA1*01:02 SUBTYPINGS BY OTHER ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No modifications made.