

101.111-24/06 – including *Taq* polymerase
101.111-24u/06u – without *Taq* polymerase

Lot No.: **19E**

Olerup SSP[®] DRB1*01

Product number: 101.111-24/06 – including *Taq* pol.
101.111-24u/06u – without *Taq* pol.
Lot number: 19E
Expiry date: 2010-March-01
Number of tests: 24 test – Product No. 101.111-24
6 tests – Product No. 101.111-06
Number of wells per test: 16

CHANGES COMPARED TO THE PREVIOUS DRB1*01 LOT (X49):

Well	5'-primer	3'-primer	rationale
7	-	Added	Primer added for the DRB1*0117 allele.
11	-	Added	Primer added for the DRB1*0118 allele.
16	New	New	New primer pair for the DRB1*0119 allele.

THE NUMBER OF WELLS has been increased from 15 to 16.

ALLELE COVERAGE:

DRB1*0101 to DRB1*0119, i.e. all the currently recognized DRB1*01 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2008-January-11, release 2.20.0.

RESOLUTION IN DRB1*01 HOMOZYGOTES:

Good.

INFLUENCE ON THE INTERPRETATION OF DRB1*01 SUBTYPINGS BY NON-DRB1*01 ALLELES:

None.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No modifications made.