

101.111-24/06 – including *Taq* polymerase  
101.111-24u/06u – without *Taq* polymerase

Lot No.: **50F**

## **Olerup SSP® DRB1\*01**

**Product number:** 101.111-24/06 – including *Taq* pol.  
101.111-24u/06u – without *Taq* pol.  
**Lot number:** 50F  
**Expiry date:** 2011-January-01  
**Number of tests:** 24 test – Product No. 101.111-24  
6 tests – Product No. 101.111-06  
**Number of wells per test:** 16

### **CHANGES COMPARED TO THE PREVIOUS DRB1\*01 LOT (19E):**

Well	5'-primer	3'-primer	rationale
9	-	Added	Primer added for the DRB1*0120 allele.

**THE NUMBER OF WELLS** is unchanged.

### **ALLELE COVERAGE:**

DRB1\*0101 to DRB1\*0120, i.e. all the currently recognized DRB1\*01 alleles, give rise to unique amplification patterns; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2009-January-16, release 2.24.0.

### **RESOLUTION IN DRB1\*01 HOMOZYGOTES:**

Good.

### **INFLUENCE ON THE INTERPRETATION OF DRB1\*01 SUBTYPINGS BY NON-DRB1\*01 ALLELES:**

None.

### **MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No modifications made.