

101.111-24/06 – including *Taq* polymerase  
101.111-24u/06u – without *Taq* polymerase

Lot No.: **28M**

## **Olerup SSP® DRB1\*01**

**Product number:** 101.111-24/06 – including *Taq* pol.  
101.111-24u/06u – without *Taq* pol.  
**Lot number:** 28M  
**Expiry date:** 2013-November-01  
**Number of tests:** 24 test – Product No. 101.111-24/24u  
6 tests – Product No. 101.111-06/06u  
**Number of wells per test:** 22

### **CHANGES COMPARED TO THE PREVIOUS DRB1\*01 LOT (90G):**

Well	5'-primer	3'-primer	rationale
20	New	New	New primer pair for the DRB1*01:33N allele.
21	New	New	New primer pairs for the DRB1*01:32 and DRB1*01:34 alleles.
22	New	New	New primer pair for the DRB1*01:36 allele.

**THE NUMBER OF WELLS** is increased from 19 to 22.

### **ALLELE COVERAGE:**

DRB1\*01:01 to DRB1\*01:36, i.e. all the currently recognized DRB1\*01 alleles, give rise to unique amplification patterns; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2011-April-08, release 3.4.0.

### **RESOLUTION IN DRB1\*01 HOMOZYGOTES:**

Good.

### **INFLUENCE ON THE INTERPRETATION OF DRB1\*01 SUBTYPINGS BY NON-DRB1\*01 ALLELES:**

None.

### **MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No modifications made.