

101.111-24/06 – including *Taq* polymerase

101.111-24u/06u – without *Taq* polymerase

Lot No.: 1S5

Olerup SSP® DRB1*01

Product number: 101.111-24/06 – including *Taq* pol.
101.111-24u/06u – without *Taq* pol.

Lot number: 1S5

Expiry date: 2027-05-01

Number of tests: 24 test – Product No. 101.111-24/24u
6 tests – Product No. 101.111-06/06u

Number of wells per test: 23+1

CHANGES COMPARED TO THE PREVIOUS DRB1*01 LOT (2R5):

The DRB1*01 primer set is unchanged compared to the previous Olerup SSP® DRB1*01 (Lot No. 2R5).

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

DRB1*01:01 to DRB1*01:141, i.e. all the currently recognized DRB1*01 alleles, will be amplified by the primers in the DRB1*01¹ subtyping kit; www.ebi.ac.uk/imgt/hla, 2023-January-12, release 3.51.0.

The DRB1*01 kit enables separation of the confirmed DRB1*01 alleles as listed in the IMGT/HLA database 3.26.0. An HLA allele is listed as confirmed by IMGT/HLA if it has been sequenced by more than a single laboratory or from multiple sources.

The DRB1*01 kit also enables identification of many null and alternatively expressed alleles.

The following DRB1*01 alleles can be distinguished by the different sizes of the HLA-specific PCR product:

Alleles	Primer mix
DRB1*01:40N, 01:67	19

¹Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

INFLUENCE ON THE INTERPRETATION OF DRB1*01 SUBTYPINGS BY NON-DRB1*01 ALLELES:

None.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.



0197

For *In Vitro* Diagnostic Use
MA125 v01 SSP Release Note Template
Date: June 2023, Rev. No: 00