

DRB1\*10  
101.129-06 – licensed for PCR  
101.129-06u – not licensed for PCR  
Lot No.: Y32

1

## Release Note - *Olerup* SSP™ DRB1\*10 – Lot No. Y32

Product number:	101.129-06 – licensed for PCR 101.129-06u – <u>not</u> licensed for PCR
Lot number:	Y32
Expiry date:	2009-October-01
Number of tests:	6
Number of tubes per test:	3

### ALLELE COVERAGE:

DRB1\*1001 to DRB1\*1002, i.e. all the currently recognized DRB1\*10 alleles, give rise to unique amplification patterns; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2007-October-05, release 2.19.0.

### RESOLUTION IN DRB1\*10 HOMOZYGOTES:

Very good.

### INFLUENCE ON THE INTERPRETATION OF DRB1\*10 SUBTYPINGS BY NON-DRB1\*10 ALLELES:

None.

### MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.