

DRB1*10
101.129-06 – including *Taq* polymerase
101.129-06u – without *Taq* polymerase
Lot No.: **84E**

Release Note

Page 1 of 1

Olerup SSP[®] DRB1*10

Product number: 101.129-06 – including *Taq* polymerase
101.129-06u – without *Taq* polymerase
Lot number: 84E
Expiry date: 2010-July-01
Number of tests: 6
Number of wells per test: 3

CHANGES COMPARED TO THE PREVIOUS DRB1*10 LOT (Y61):

The DRB1*10 primer set is unchanged.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

DRB1*1001 to DRB1*1002, i.e. all the currently recognized DRB1*10 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2008-April-08, release 2.21.0.

RESOLUTION IN DRB1*10 HOMOZYGOTES:

Very good.

INFLUENCE ON THE INTERPRETATION OF DRB1*10 SUBTYPINGS BY NON-DRB1*10 ALLELES:

None.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.