

## Olerup SSP<sup>®</sup> DRB1\*10

**Product number:** 101.129-06 – including *Taq* polymerase  
101.129-06u – without *Taq* polymerase  
**Lot number:** 81F  
**Expiry date:** 2011-April-01  
**Number of tests:** 6  
**Number of wells per test:** 4

### CHANGES COMPARED TO THE PREVIOUS DRB1\*10 LOT (84E):

Well	5'-primer	3'-primer	rationale
4	New	New	New primer pair for the DRB1*1003 allele.

**THE NUMBER OF WELLS** has been increased from 3 to 4.

### ALLELE COVERAGE:

DRB1\*1001 to DRB1\*1003, i.e. all the currently recognized DRB1\*10 alleles, give rise to unique amplification patterns; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2009-April-17, release 2.25.0.

### RESOLUTION IN DRB1\*10 HOMOZYGOTES:

Very good.

### INFLUENCE ON THE INTERPRETATION OF DRB1\*10 SUBTYPINGS BY NON-DRB1\*10 ALLELES:

None.

### MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.