

DRB1*10
101.129-06 – including *Taq* polymerase
101.129-06u – without *Taq* polymerase
Lot No.: **79G**

Release Note

Page 1 of 1

Olerup SSP[®] DRB1*10

Product number: 101.129-06 – including *Taq* polymerase
101.129-06u – without *Taq* polymerase
Lot number: 79G
Expiry date: 2012-January-01
Number of tests: 6
Number of wells per test: 4

CHANGES COMPARED TO THE PREVIOUS DRB1*10 LOT (81F)

The DRB1*10 primer set is unchanged compared to the previous lot.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

DRB1*1001 to DRB1*1003, i.e. all the currently recognized DRB1*10 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2009-October-19, release 2.27.0.

RESOLUTION IN DRB1*10 HOMOZYGOTES:

Very good.

INFLUENCE ON THE INTERPRETATION OF DRB1*10 SUBTYPINGS BY NON-DRB1*10 ALLELES:

None.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.