

101.129-06 – including *Taq* polymerase101.129-06u – without *Taq* polymerase

Lot No.: 13L

Olerup SSP® DRB1*10

Product number: 101.129-06 – including *Taq* polymerase
101.129-06u – without *Taq* polymerase
Lot number: 13L
Expiry date: 2013-August-01
Number of tests: 6
Number of wells per test: 4

CHANGES COMPARED TO THE PREVIOUS DRB1*10 Lot (79G)

The DRB1*10 primer set is unchanged compared to the previous lot.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

DRB1*10:01 to DRB1*10:03, i.e. all the currently recognized DRB1*10 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2010-October-15, release 3.2.0.

RESOLUTION IN DRB1*10 HOMOZYGOTES:

Very good.

INFLUENCE ON THE INTERPRETATION OF DRB1*10 SUBTYPINGS BY NON-DRB1*10 ALLELES:

None.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.