

Olerup SSP[®] DRB1*12

Product number: 101.128-12 – including *Taq* polymerase
101.128-12u – without *Taq* polymerase
Lot number: 25E
Expiry date: 2010-March-01
Number of tests: 12
Number of wells per test: 22

CHANGES COMPARED TO THE PREVIOUS DRB1*12 LOT (X21):

Well	5'-primer	3'-primer	rationale
22	New	New	New primer pair for the DRB1*1216 allele.

THE NUMBER OF WELLS has been increased from 21 to 22.

ALLELE COVERAGE:

DRB1*1201 to DRB1*1216, i.e. all the currently recognized DRB1*12 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2008-January-11, release 2.20.0.

RESOLUTION IN DRB1*12 HOMOZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF DRB1*12 SUBTYPINGS BY NON-DRB1*12 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.