

## **Olerup SSP<sup>®</sup> DRB1\*12 Add-on**

<b>Product number:</b>	<b>101.815-12 – including <i>Taq</i> polymerase 101.815-12u – without <i>Taq</i> polymerase</b>
<b>Lot number:</b>	<b>51S</b>
<b>Expiry date:</b>	<b>2015-December-01</b>
<b>Number of tests:</b>	<b>12</b>
<b>Number of wells per test:</b>	<b>2</b>

### **CHANGES COMPARED TO THE PREVIOUS DRB1\*12 ADD-ON LOT (31N):**

The DRB1\*12 Add-on primer set is unchanged compared to the previous lot.

**THE NUMBER OF WELLS** is 2 in an 8 well plate.

### **ALLELE COVERAGE:**

The DRB1\*12:01 and DRB1\*12:10 alleles give different patterns in the DRB1\*12 Add-on subtyping kit<sup>1</sup>; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2013-April-17, release 3.12.0.

<sup>1</sup>The DRB1\*12 Add-on subtyping kit cannot distinguish the following silent mutations: the DRB1\*12:01:01-12:01:05 alleles.

### **RESOLUTION IN HOMO- AND HETEROZYGOTES:**

Excellent.

### **INFLUENCE ON THE INTERPRETATION OF DRB1\*12 ADD-ON SUBTYPINGS BY NON-DRB1\*12 ALLELES:**

None.