

101.815-12 – including *Taq* polymerase

101.815-12u – without *Taq* polymerase

Lot No.: **9F3**

Olerup SSP[®] DRB1*12 Add-on

Product number:	101.815-12 – including <i>Taq</i> polymerase 101.815-12u – without <i>Taq</i> polymerase
Lot number:	9F3
Expiry date:	2020-06-01
Number of tests:	12
Number of wells per test:	3+1

The DRB1*12 Add-on primer set is unchanged compared to the previous Olerup SSP[®] DRB1*12 Add-on (Lot No. 9D0).

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

The DRB1*12:01, DRB1*12:10 and DRB1*12:17 alleles give different patterns in the DRB1*12 Add-on subtyping kit¹; www.ebi.ac.uk/imgt/hla, 2017-August-10, release 3.29.0.

¹Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

RESOLUTION IN HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF DRB1*12 ADD-ON SUBTYPINGS BY NON-DRB1*12 ALLELES:

None.