

101.815-12 – including *Taq* polymerase

101.815-12u – without *Taq* polymerase

Lot No.: 7G3

Olerup SSP® DRB1*12 Add-on

Product number:	101.815-12 – including <i>Taq</i> polymerase 101.815-12u – without <i>Taq</i> polymerase
Lot number:	7G3
Expiry date:	2021-01-01
Number of tests:	12
Number of wells per test:	3+1

The DRB1*12 Add-on primer set is unchanged compared to the previous *Olerup* SSP® DRB1*12 Add-on (**Lot No. 9F3**).

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

The DRB1*12:01, DRB1*12:10 and DRB1*12:17 alleles give different patterns in the DRB1*12 Add-on subtyping kit¹; www.ebi.ac.uk/imgt/hla, 2018-April-16, release 3.32.0.

¹Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

RESOLUTION IN HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF DRB1*12 ADD-ON SUBTYPINGS BY NON-DRB1*12 ALLELES:

None.