

101.125-24/06 – including *Taq* polymerase
101.125-24u/06u – without *Taq* polymerase

Lot No.: **55F**

Olerup SSP® DRB1*15

Product number: 101.125-24/06 – including *Taq* pol.
101.125-24u/06u – without *Taq* pol.
Lot number: 55F
Expiry date: 2011-February-01
Number of tests: 24 test – Product No. 101.125-24
6 tests – Product No. 101.125-06
Number of wells per test: 24

CHANGES COMPARED TO THE PREVIOUS DRB1*15 LOT (35E):

Well	5'-primer	3'-primer	rationale
18	-	-	Changed positive control primer pair
19	-	Added	New primer for the DRB1*1530 allele.
20	-	Added	New primer for the DRB1*1527 allele, for improved resolution of the DRB1*1527 and DRB1*1531 alleles.
22	Added	-	New primer for the DRB1*1528 allele.
23	-	Added	New primer for the DRB1*1529 allele.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

DRB1*1501 to DRB1*1531, i.e. all the currently recognized DRB1*15 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2000-January-16, release 2.24.0.

RESOLUTION IN DRB1*15 HOMOZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF DRB1*15 SUBTYPINGS BY NON-DRB1*15 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.