

DRB1\*15/16 low resolution  
101.151-24 – licensed for PCR  
101.151-24u – not licensed for PCR  
Lot No.: Y18

1

## **Olerup SSP™ DRB1\*15/16 low resolution**

Product number: 101.151-24 – licensed for PCR  
101.151-24u – not licensed for PCR  
Lot number: Y18  
Expiry date: 2009-July-01  
Number of tests: 24  
Number of tubes per test: 2

### **CHANGES COMPARED TO THE PREVIOUS DRB1\*15/16 LOW RESOLUTION LOT (R95):**

Tube	5'-primer	3'-primer	rationale
1	-	Exchanged	Primer exchanged to amplify all DRB1*15 alleles and no DRB1*16 alleles.
2	-	Exchanged	Primer exchanged to amplify all DRB1*16 alleles and no DRB1*15 alleles.

**THE NUMBER OF WELLS** is unchanged.

### **ALLELE COVERAGE:**

DRB1\*1501 to DRB1\*1522 and DRB1\*1601 to DRB1\*1611, i.e. all the currently recognized DRB1\*15 and DRB1\*16 alleles will be amplified by the primers in the DRB1\*15/16 subtyping kit and grouped into the DRB1\*15 and DRB1\*16 groups, respectively; 2007-July-09, release 2.18.0.

### **INFLUENCE ON THE INTERPRETATION OF DRB1\*15/16 LOW RESOLUTION SUBTYPINGS BY NON-DRB1\*15/16 ALLELES:**

None.

### **MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No modifications made.