

Lot No.: **48E**

## **Olerup SSP<sup>®</sup> DRB1\*15/16 low resolution**

**Product number:** 101.151-24 – including *Taq* polymerase  
101.151-24u – without *Taq* polymerase  
**Lot number:** 48E  
**Expiry date:** 2010-May-01  
**Number of tests:** 24  
**Number of wells per test:** 2

### **CHANGES COMPARED TO THE PREVIOUS DRB1\*15/16 LOW RESOLUTION LOT (Y18):**

<b>Well</b>	<b>5'-primer</b>	<b>3'-primer</b>	<b>rationale</b>
1	-	Added	Primer added for the DRB1*1527 allele.
2	Added	-	Primer added for the DRB1*1612 allele.

**THE NUMBER OF WELLS** is unchanged.

#### **ALLELE COVERAGE:**

DRB1\*1501 to DRB1\*1527 and DRB1\*1601 to DRB1\*1613N, i.e. all the currently recognized DRB1\*15 and DRB1\*16 alleles will be amplified by the primers in the DRB1\*15/16 subtyping kit and grouped into the DRB1\*15 and DRB1\*16 groups, respectively; 2008-April-08, release 2.21.0.

#### **INFLUENCE ON THE INTERPRETATION OF DRB1\*15/16 LOW RESOLUTION SUBTYPINGS BY NON-DRB1\*15/16 ALLELES:**

None.

#### **MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No modifications made.