

101.126-12 – including *Taq* polymerase101.126-12u – without *Taq* polymeraseLot No.: **02M****Olerup SSP[®] DRB1*16**

Product number: 101.126-12 – including *Taq* polymerase
101.126-12u – without *Taq* polymerase
Lot number: 02M
Expiry date: 2013-October-01
Number of tests: 12
Number of wells per test: 16

CHANGES COMPARED TO THE PREVIOUS DRB1*16 LOT (04K):

Well	5'-primer	3'-primer	rationale
12	Added	-	Primer added for the DRB1*16:17 allele.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

DRB1*16:01 to DRB1*16:18, i.e. all the currently recognized DRB1*16 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2011-January-14, release 3.3.0.

RESOLUTION IN DRB1*16 HOMOZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF DRB1*16 SUBTYPINGS BY NON-DRB1*16 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.