

101.113-24/04 – including *Taq* polymerase  
101.113-24u/04u – without *Taq* polymerase

Lot No.: **52M**

## **Olerup SSP® DRB1\*03**

Product number:	101.113-24/04 – including <i>Taq</i> pol. 101.113-24u/04u – without <i>Taq</i> pol.
Lot number:	52M
Expiry date:	2014-January-01
Number of tests:	24 tests – Product No. 101.113-24 4 tests – Product No. 101.113-04
Number of wells per test:	32

### **CHANGES COMPARED TO THE PREVIOUS DRB1\*03 LOT (01K):**

The DRB1\*03 primer set is unchanged compared to the previous lot.

**THE NUMBER OF WELLS** is unchanged.

### **ALLELE COVERAGE:**

DRB1\*03:01 to DRB1\*03:52, i.e. all the currently recognized DRB1\*03 alleles, give rise to unique amplification patterns; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2010-April-01, release 3.0.0.

### **RESOLUTION IN DRB1\*03 HOMOZYGOTES:**

The DRB1\*03:01,03:01, DRB1\*03:01,03:02 and DRB1\*03:02,03:02 genotypes all give rise to unique amplification patterns.

### **INFLUENCE ON THE INTERPRETATION OF DRB1\*03 SUBTYPINGS BY NON-DRB1\*03 ALLELES:**

Most frequently encountered DRB1\*03,11; DRB1\*03,13 and DRB1\*03,14 genotypes give rise to unique amplification patterns.

### **MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No comments received.