

Lot No.: **40M**

Olerup SSP[®] DRB1*04

Product number: 101.114-24/03 – including *Taq* pol.
101.114-24u/03u – without *Taq* pol.
Lot number: 40M
Expiry date: 2014-February-01
Number of tests: 24 tests – Product No. 101.114-24
3 tests – Product No. 101.114-03
Number of wells per test: 45

CHANGES COMPARED TO THE PREVIOUS DRB1*04 LOT (28K):

Well	5'-primer	3'-primer	rationale
23	Moved	Moved	One primer pair moved to well 45.
45	New	New	Primer pair from well 23.

THE NUMBER OF WELLS has been increased from 44 to 45

ALLELE COVERAGE:

DRB1*04:01 to DRB1*04:89, i.e. all the currently recognized DRB1*04 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2010-April-01 release 3.0.0.

RESOLUTION IN DRB1*04 HOMOZYGOTES:

All possible homo- and heterozygous genotypic combinations between the DRB1*04:01 to 04:10 alleles give rise to unique amplification patterns.

INFLUENCE ON THE INTERPRETATION OF DRB1*04 SUBTYPINGS BY NON-DRB1*04 ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.