

101.814-12 – including *Taq* polymerase

101.814-12u – without *Taq* polymerase

Lot No.: **4F1**

## Olerup SSP<sup>®</sup> DRB1\*04 Add-on

Product number:	101.814-12 – including <i>Taq</i> polymerase 101.814-12u – without <i>Taq</i> polymerase
Lot number:	4F1
Expiry date:	2019-12-01
Number of tests:	12
Number of wells per test:	3+1

The DRB1\*04 Add-on primer set is unchanged compared to the previous *Olerup* SSP<sup>®</sup> DRB1\*04 Add-on (**Lot No. 7D0**).

**THE NUMBER OF WELLS** is unchanged.

### ALLELE COVERAGE:

The DRB1\*04:07 and DRB1\*04:92 alleles give different patterns in the DRB1\*04 Add-on subtyping kit<sup>1</sup>; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2016-October-14, release 3.26.0.

The primer set also resolves the SBT heterozygous ambiguities:

DRB1\*04:02:01, 04:03:01 = DRB1\*04:37, 04:88

DRB1\*04:02:01, 04:04:01 = DRB1\*04:37, 04:56

<sup>1</sup>Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

### RESOLUTION IN HOMO- AND HETEROZYGOTES:

Excellent.

### INFLUENCE ON THE INTERPRETATION OF DRB1\*04 ADD-ON SUBTYPINGS BY NON-DRB1\*04 ALLELES:

None.